

**Esthetic Restoration Matrix**  
**Rocky Mountain Dental Laboratory, Inc.**

Category	Indications	Advantages	Disadvantages
<i>Feldspathic Ceramics</i>	Veneers All-ceramic crowns	<ul style="list-style-type: none"> <li>• Natural translucency &amp; vitality</li> <li>• Can be fabricated as thin as 0.2 mm</li> <li>• Can be designed almost clear if the natural underlying tooth is close to final desired shade.</li> <li>• Ability to alter translucencies &amp; opacities within single restoration to block out discolored areas &amp; to allow ideal tooth prep to shine through restoration.</li> <li>• Marginal fit &amp; accuracy is excellent, especially @ gingival margins.</li> </ul>	<ul style="list-style-type: none"> <li>▪ Must be adhesively bonded to achieve adequate strength.</li> <li>▪ Uniform ceramic thickness causes weakness when the thickness is greater than 2.5 mm. Can require buildup, core or base to provide most ideal preparation</li> </ul>
<i>Lithium Disilicate (e.max CAD)</i>	Anterior/posterior restorations, Veneers, partial crowns	<ul style="list-style-type: none"> <li>• Any cementation technique is acceptable.</li> <li>• Stronger than first generation of millable feldspar ceramics (360μPa vs. 175μPa).</li> <li>• Very enamel friendly &amp; great esthetics.</li> </ul>	<ul style="list-style-type: none"> <li>▪ Contraindications are full veneers on molar crowns, deep sub-gingival preparations, patients with severely reduced residual dentitions and bruxism. Flexural strength is lower (360μPa).</li> </ul>
<i>Yttrium stabilized-Zirconia (YZ)</i>	Bridges up to 7 units All-ceramic restorations	<ul style="list-style-type: none"> <li>• Very dense, strong framework (960μPa)</li> <li>• Final restorations are consistently higher than metal supported restoration</li> <li>• Conventional cement. No Etching.</li> </ul>	<ul style="list-style-type: none"> <li>▪ NONE AT THIS TIME</li> </ul>